

GGIR package

Study aim and background

- Sickness behaviors, including fatigue, depression, and disruptions in sleep and activity, are among the most common and distressing side effects of cancer and its treatment.
- To date, research has focused on identifying the prevalence and severity of sickness behaviors at various points along the cancer continuum and establishing the efficacy of behavioral interventions to reduce sickness behaviors.
- The next generation of research should refine our current understanding of how sickness behaviors develop.

Study aim and background

- Previous research to examine onset of sickness behaviors in women undergoing platinum-based chemotherapy for gynecologic cancer, one of the most arduous treatment regimens for cancer.
 - Found that sickness behaviors previously thought to cluster together actually occur in a cascade pattern
 - disruptions in sleep contributing to lagged increases in fatigue and fatigue in turn contributing to lagged increases in depression.
- This work is clinically important because it suggests that interventions targeted to sickness behaviors early in the cascade may have beneficial effects on multiple downstream sickness behaviors.


Study aim and background

- Recruit a sample of 150 women receiving chemotherapy for gynecologic cancer and a matched comparison sample of 150 women with no history of cancer.
- Both samples will be evaluated on eight occasions during which sickness behaviors (and pro-inflammatory cytokines) will be assessed.
- Assessments will occur before, during, and at six and twelve months after chemotherapy.

Study aim and background

- **To identify the magnitude and features of changes in sickness behaviors (i.e., fatigue, depression, and disruptions in sleep and activity) and pro-inflammatory cytokines in gynecologic cancer patients treated with chemotherapy compared to women without cancer.**
 - The cancer sample is expected to demonstrate greater increases in sickness behaviors during chemotherapy than the non-cancer sample.
 - Comparisons will also provide information on whether sickness behaviors in the cancer sample are elevated prior to chemotherapy and whether they increase, decrease, or remain stable following completion of chemotherapy.
 - Sickness behaviors in patients are expected to increase in a cascade pattern, with sleep disruptions predicting lagged increases in fatigue and fatigue predicting lagged increases in depression.

Timing of Assessments



GGIR

- A tool to process and analyze data collected with wearable raw acceleration sensors [as described in van Hees and colleagues (2014)]
- Developed and tested for binary data from 'GENEActiv' and GENEActiv devices, .csv-export data from 'Actigraph' devices, and .cwa and .wav-format data from 'Axivity'

----- Data File Created By ActiGraph GT3X+ ActiLife v6.13.3 Firmware v2.2.1 date format M/d/yyyy at
30 Hz Filter Normal -----
Serial Number: MRA3B09130251
Start Time 14:00:00
Start Date 5/24/2013
Epoch Period (hh:mm:ss) 00:00:00
Download Time 12:48:47
Download Date 6/7/2013
Current Memory Address: 0
Current Battery Voltage: 3.99 Mode = 12

Accelerometer X,Accelerometer Y,Accelerometer Z

1.079,-0.176,-0.126

1.05,-0.155,-0.076

1.012,-0.155,-0.029

1.009,-0.15,0

0.997,-0.158,0.044

1.041,-0.147,0.135

1.097,-0.129,0.176

1.047,-0.12,0.185

1.018,-0.114,0.17

1.032,-0.114,0.17

1.062,-0.114,0.147

0.982,-0.085,0.15

0.965,-0.073,0.176

.

.

.

.

14,000,000 plus.... rows

Main GGIR functions

The function arguments need to be tailored to your experimental protocol.

GGIR is structured in 5 parts.

- Part 1: Loads the data and stores derived features (aggregations) needed for the other parts. This is the time consuming part.
- Part 2: Data quality analyses and low level description of signal features per person per day and per file. At this point a day is defined from midnight to midnight
- Part 3: Estimation of sustained inactivity and sleep periods, needed for input to Part 4 for sleep detection
- Part 4: Labels the sustained inactive periods detected in Part 3 as sleep, or daytime sustained inactivity, per night and per file
- Part 5: Derives sleep and physical activity characteristics by re-using information derived in part 2, 3 and 4. Total time in intensity categories, the number of bouts, time spent in bouts and average acceleration (overall activity) is calculated.

Part 1

- Part 1: Loads the data, creates file structure, and stores derived features (aggregations) needed for the other parts. This is the time consuming part. Once this is done, parts 2-5 can be run (or re-run with different parameters in parts 2-5) relatively quickly.

Part 1

- Calls function `g.getmeta` and `g.calibrate`, and converts the output to `.RData`-format which will be the input for `g.part2`
- Function generates a folder structure to keep track of various output files.
- `g.part 1` takes much longer to run than other parts and involves only minor decisions of interest to the movement scientist..

codetalker - zaza - SSH Secure File Transfer

File Edit View Operation Window Help

Quick Connect Profiles

Remote Folders

- Output
 - OutputWITHlog010719
 - OutputWITHlog052019
 - output_Period1
 - meta
 - basic
 - ms2.out
 - ms3.out
 - ms4.out
 - ms5.out
 - sleep.qc
 - results
 - QC
 - output_Period2
 - meta
 - results
 - output_Period3
 - meta
 - results
 - output_Period4
 - meta
 - results
 - output_Period5
 - output_Period6
 - output_Period7
 - output_Period8

Remote Name	Size	Type	Modified	Attribu...
basic		Folder	05/21/2019 11:...	drwxr-...
ms2.out		Folder	05/21/2019 04:...	drwxr-...
ms3.out		Folder	05/21/2019 06:...	drwxr-...
ms4.out		Folder	05/22/2019 07:...	drwxr-...
ms5.out		Folder	05/22/2019 08:...	drwxr-...
sleep.qc		Folder	05/21/2019 06:...	drwxr-...
sleeplog.RData	4,443	RDATA...	05/22/2019 07:...	-rwxr-x...

Transfer Queue

codetalker - zaza - SSH Secure File Transfer

File Edit View Operation Window Help

Quick Connect Profiles

Remote Folders

- Output
 - OutputWITHlog010719
 - OutputWITHlog052019
 - output_Period1
 - meta
 - basic
 - ms2.out
 - ms3.out
 - ms4.out
 - ms5.out
 - sleep.qc
 - results
 - QC
 - output_Period2
 - meta
 - results
 - output_Period3
 - meta
 - results
 - output_Period4
 - meta
 - results
 - output_Period5
 - output_Period6
 - output_Period7
 - output_Period8

Remote Name / Size Type Modified Attribu...

QC		Folder	05/22/2019 08:...	drwxr-...
part2_daysummary.csv	1,224,5...	Micros...	05/22/2019 08:...	-rwxr-x...
part2_summary.csv	962,650	Micros...	05/22/2019 08:...	-rwxr-x...
part4_nightsunmary_sleep_cleaned.csv	202,404	Micros...	05/22/2019 08:...	-rwxr-x...
part4_summary_sleep_cleaned.csv	146,733	Micros...	05/22/2019 08:...	-rwxr-x...
part5_daysummary_WW_L20M60V240_T5A5.csv	1,113,8...	Micros...	05/22/2019 08:...	-rwxr-x...
part5_personsummary_WW_L20M60V240_T5A5.csv	353,707	Micros...	05/22/2019 08:...	-rwxr-x...
visualisation_sleep.pdf	396,747	Adobe ...	05/22/2019 07:...	-rwxr-x...

codetalker - zaza - SSH Secure File Transfer

File Edit View Operation Window Help

Quick Connect Profiles

Remote Folders

- Output
 - OutputWITHlog010719
 - OutputWITHlog052019
 - output_Period1
 - meta
 - basic
 - ms2.out
 - ms3.out
 - ms4.out
 - ms5.out
 - sleep.qc
 - results
 - QC
 - output_Period2
 - meta
 - results
 - output_Period3
 - meta
 - results
 - output_Period4
 - meta
 - results
 - output_Period5
 - output_Period6
 - output_Period7
 - output_Period8

Remote Name / Size Type Modified Attribu...

data_quality_report.csv	37,192	Micros...	05/22/2019 08:...	-rwxr-x...
filelocationkey.csv	1,290	Micros...	05/22/2019 09:...	-rwxr-x...
part4_nightssummary_sleep_full.csv	246,081	Micros...	05/22/2019 08:...	-rwxr-x...
part4_summary_sleep_full.csv	159,140	Micros...	05/22/2019 08:...	-rwxr-x...
plots_to_check_data_quality_1.pdf	880,146	Adobe ...	05/22/2019 08:...	-rwxr-x...
sessioninfo_part1.RData	6,051	RDATA...	05/22/2019 09:...	-rwxr-x...
sessioninfo_part2.RData	6,055	RDATA...	05/22/2019 08:...	-rwxr-x...
sessioninfo_part3.RData	6,055	RDATA...	05/22/2019 07:...	-rwxr-x...
sessioninfo_part4.RData	6,055	RDATA...	05/22/2019 08:...	-rwxr-x...

Part 1

- `g.getmeta` Reads an accelerometer file in blocks, extracts various features and stores average feature value per short or long epoch.
 - Acceleration and angle metrics are stored at short epoch length.
 - The non-wear indication score, the clipping score, temperature (if available), light (if available), and Euclidean norm are stored at long epoch length.
- `g.calibrate` Function starts by identifying ten second windows of non-movement. Next, the average acceleration per axis per window is used to estimate calibration error (offset and scaling) per axis. The function provides recommended correction factors to address the calibration error and a summary of the calibration procedure.

g.Part1 key parameters

- **datadir** = Directory where the accelerometer files are stored or list of accelerometer file-names and directories
- **Outputdir** = where the output needs to be stored. Note that this function will attempt to create folders in this directory and uses those folder to organise output
- **f0**=1, File index to start with (default = 1). Index refers to the filenames sorted in increasing order
- **F1** f1File index to finish with (defaults to number of files availa

g.Part1 key parameters

- **do.bfen** if TRUE, calculate metric BFEN with band-pass filter configuration set by lb and hb
- **do.enmo** if TRUE (default), calculate metric ENMO
- **do.lfenmo** if TRUE, calculate metric LFENMO with low-pass filter
 - **do.hfen**
 - **do.hfenplus**
- **do.mad** if TRUE, calculate metric MAD (Mean Amplitude Deviation)
- **mvpthreshold** Threshold for MVPA estimation. Threshold needs to be based on metric ENMO. This can be a single number or an vector of numbers, e.g. c(100,120). In the latter case the code will estimate MVPA separately for each threshold.
- **do.cal** Whether to apply auto-calibration or not, see g.calibrate. Default and recommended setting is TRUE

g.Getmeta key parameters

- **windowsizes**

- Three values to indicate the lengths of the windows as in `c(window1,window2,window3)`:
- `window1` is the short epoch length in seconds and by default 5 this is the time window over which acceleration and angle metrics are calculated,
- `window2` is the long epoch length in seconds for which non-wear and signal clipping are defined, default 900.
- `window3` is the window length of data used for non-wear detection and by default 3600 seconds.

Metrics

- **Euclidean Norm Minus One**
- **Mean Amplitude Deviation**
- **Others:** BFEN with band-pass filter
 HFEN with low-pass filter
 etc... etc...
- **No one metric systematically outperforms all others across a wide range of standardized kinematic conditions.**
- **ENMO and MAD and certain thresholds in Bakrania et al. 2016 are robust for separating sedentary behaviors from light activity**

ENMO

$$r_i = \sqrt{x_i^2 + y_i^2 + z_i^2} - 1000$$

1000 = 1000 milligravitational units = 1 g

MAD

$$MAD = \frac{1}{n} \sum_{i=1}^n |r_i - \bar{r}|$$

Where $r_i = \sqrt{x_i^2 + y_i^2 + z_i^2}$

Is the i^{th} vector magnitude at each time point

and $\bar{r} =$

mean vector magnitude, $n = \text{length of time}$

g.Calibrate key parameters

- **chunksize** number between 0.2 and 1 to specify the size of chunks to be loaded as a fraction of a 12 hour period, e.g. 0.5 equals 6 hour chunks. The default is 1 (12 hrs). For machines with less than 4Gb of RAM memory a value below 1 is recommended.
- **minloadcrit** the minimum number of hours the code needs to read for the autocalibration procedure to be effective. After loading these hours only extra data is loaded if calibration error has not been reduced to under 0.01 g.





Autocalibration

An acceleration sensor works on the principle that acceleration is captured mechanically and converted into an electrical signal.

The relationship between the electrical signal and the acceleration is usually assumed to be linear, involving an offset and a gain factor.

Autocalibration

- A recording of acceleration is screened for non-movement periods.
- Next, the moving average over the non-movement periods is taken from each of the three orthogonal sensor axes and used to generate a three-dimensional ellipsoid representation that should ideally be a sphere with radius 1 g
 - g- units (1 standard g = 9.80665 m/s^2)
- Deviations between the radius of the three-dimensional ellipsoid and 1 g (ideal calibration) can then be used to derive correction factors for sensor axis-specific calibration error





Part 2

- Loads the output from g.part1 and then applies g.impute and g.analyse, after which the output is converted to .RData-format which will be used by g.shell.GGIR to generate reports.
- Part 2 generates:
 - part2_summary.csv
 - part2_daysummary.csv
 - QC/data_quality_report.csv
 - QC/plots to check data quality 1.pdf


g.Impute key parameters


- **strategy** how to deal with knowledge about study protocol.
 - value = 1 to select data based on hrs.del.start, hrs.del.end, and maxdur.
 - Value = 2 to only use the data between the first midnight and the last midnight
 - value = 3 only selects the most active X days in the files. X is specified by argument ndayswindow
- Others depend on value of **strategy**






device: actigraph | filename: 0108-1.csv







Part 3

- Estimates the sustained inactivity periods in each day
- Used as input for g.part4 which then labels them as nocturnal sleep or day time sustained inactivity periods.

Part 3 key parameters

- `anglethreshold` Angle threshold (degrees) for sustained inactivity periods detection, default = 5
- `timethreshold` Time threshold (minutes) for sustained inactivity periods detection, default = 5. This can be specified as multiple thresholds, each of which will be implemented. For example, `timethreshold = c(5,10)`
- `ignorenonwear` If TRUE then ignore detected monitor non-wear periods to avoid confusion between monitor non-wear time and sustained inactivity (default = TRUE)

Part 4

- Loads output from g.part3 as stored in milestone data and sleep log information (if available)
- Uses these information sources to define nocturnal sleep and daytime sustained inactivity
- Part 4 generates the following output:
 - part4_nightsummary_sleep_cleaned.csv
 - part4_summary_sleep_cleaned.csv
 - visualisation_sleep.pdf (optional)


Part 4

- **Loglocation** Location of the spreadsheet (csv) with sleep log information.
- **outliers.only** = FALSE will visualize all subjects, if TRUE will visualize only for nights with a difference in onset or waking time larger than the variable of argument **criterror**.
- **relyonsleeplog** = If TRUE then sleep onset and waking time are defined based on timestamps derived from sleep log if FALSE (default) the sleep log is only used to guide the accelerometer-based detection. If participants were instructed NOT to wear the accelerometer during waking hours then set to TRUE, in all other scenarios set to FALSE.
- **do.visual** = TRUE generate a pdf with a visual representation of the overlap between the sleeplog entries and the accelerometer detections.
- **nnights** = Number of nights for which sleep log information should be available. It assumes that this is constant within a study. If sleep log information is missing for certain nights then leave these blank

Sleep Log

- The spreadsheet needs to have the following structure: one column for participant id, and then followed by alternatingly one column for onset time and one column for waking time (see example below). There can be multiple sleeplogs in the same spreadsheet. The first row of the spreadsheet needs to be filled with column names, it does not matter what these column names are. Timestamps are to be stored without date as in hh:mm:ss. If onset corresponds to lights out or intention to fall asleep, then it is the end-users responsibility to account for this in the interpretation of the results.

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	
1	STNO	FGASLP2	FGAWK2	FGASLP3	FGAWK3	FGASLP4	FGAWK4	FGASLP5	FGAWK5	FGASLP6	FGAWK6	FGASLP7	FGAWK7	FGASLP8	FGAWK8	FGASLP9	FGAWK9	FGASLP10	FGAWK10	
2	10002	21:40:00	5:35:00	23:30:00	6:45:00	22:20:00	6:30:00	22:45:00	7:45:00	22:30:00	6:20:00	22:20:00	6:30:00	22:20:00	6:10:00	22:30:00	6:40:00	23:25:00	6:55:00	
3	10004	23:30:00	7:00:00	22:50:00	7:20:00	22:45:00	7:45:00	22:30:00	7:30:00	22:15:00	6:45:00	22:45:00	6:45:00	22:45:00	6:45:00	23:00:00	6:30:00	23:00:00	7:00:00	
4	10005	22:15:00	6:45:00	23:00:00	6:45:00	22:30:00	8:00:00	20:00:00	7:00:00	21:00:00	6:30:00	22:00:00	6:45:00	22:00:00	6:45:00	22:00:00	7:15:00	22:00:00	7:30:00	
5	10011	0:00:00	7:25:00	23:30:00	7:25:00	23:30:00	7:25:00	23:15:00	7:30:00	23:15:00	7:30:00	23:30:00	7:45:00	0:00:00	8:30:00	23:35:00	7:30:00	23:15:00	7:35:00	
6	10013	23:00:00	6:23:00	22:50:00	6:30:00	22:50:00	6:45:00	22:50:00	6:20:00	23:00:00	6:40:00	22:50:00	6:30:00	22:55:00	6:30:00	23:00:00	6:30:00	23:10:00	6:40:00	
7	10016	23:30:00	7:10:00	0:00:00	6:30:00	23:45:00	7:00:00	23:45:00	6:00:00	23:30:00	6:45:00	0:00:00	6:15:00	0:15:00	6:50:00	23:50:00	5:50:00	23:45:00	6:15:00	
8	10017	23:00:00	8:00:00	23:00:00	8:15:00	22:30:00	9:00:00	0:00:00	9:00:00	20:22:00	8:00:00	22:30:00	8:15:00	22:15:00	8:15:00	23:30:00	8:00:00	23:00:00	8:30:00	



Some output variables

- **acc_onset** Detected onset of sleep expressed as hours since the midnight of the previous night.
- **acc_wake** Detected waking time (after sleep period) expressed as hours since the midnight of the previous night

Some output variables

- `acc_SptDuration` = Difference between onset and waking time.
- `acc_def` = Definition of sustained inactivity by accelerometer
- `sleeplog_onset` Start of Sleep Period Time window derived from (in order of priority) `sleeplog`, detected by the HDCZA algorithm, detected by L5HR6 algorithm, or specified by researcher,
- `sleeplog_wake` End of Sleep Period Time window derived from (in order of priority) `sleeplog`, detected by the HDCZA algorithm, detected by L5HR6 algorithm, or specified by resear

Part 5

- Merges the output from g.part2 and g.part4 into one report enhanced with profiling of sleep and physical activity stratified across intensity levels

Part 5 key parameters

- **threshold.lig** = Threshold for light physical activity to separate inactivity from light activity. Value can be one number or an vector of multiple numbers, e.g. `threshold.lig = c(30,40)`. If multiple numbers are entered then analysis will be replicated for each combination of threshold values.
 - **threshold.mod** (similar)
 - **threshold.vig** (similar)
- Bouts of MVPA are identified as all 5- or 10-min time windows that start with a 5-s epoch value equal to or higher than ie. 100 mg and for which 80% of subsequent 5-s epoch values are equal to or higher than the 100-mg threshold
- **boutcriter.mvpa** = A number between 0 and 1 and defines what fraction of a bout needs to be above the mvpa threshold
 - **boutcriter.in**
 - **boutcriter.lig**

Part 5 key parameters


- **boutcriter.mvpa** = Durations of mvpa bouts in minutes to be extracted. The default values are c(1,5,10) and will start with the identification of 10 minute bouts, followed by 5 minute bouts in the rest of the data, and followed by 1 minute bouts in the rest of the data.
- **boutdur.in** = Durations of inactivity bouts in minutes to be extracted. Inactivity bouts are detected in the segments of the data which were not labelled as sleep or MVPA bouts.
 - **boutdur.lig**
 - **boutdur.mvpa**
- **timewindow** = Time window over which summary statistics are derived. Value can be "MM"(midnight to midnight), "WW" (waking time to waking time)

Day 1: Wednesday 6/3/2013

What we think you did:

Your arm angle (up or down):

Your arm movement:




Day 2: Thursday 7/3/2013

What we think you did:

Your arm angle (up or down):

Your arm movement:




Day 3: Friday 8/3/2013

What we think you did:

Your arm angle (up or down):

Your arm movement:




Day 4: Saturday 9/3/2013

What we think you did:

Your arm angle (up or down):

Your arm movement:




Day 5: Sunday 10/3/2013

What we think you did:

Your arm angle (up or down):

Your arm movement:



- <https://cran.r-project.org/web/packages/GGIR/vignettes/GGIR.html#introduction>

```
mainscript_GGIR_demo.R
Source on Save Run Source
9 outputdir="C:/Users/FlorianHuber/OneDrive - Netherlands eScience Center/Sprint_project_GGIR/Test-Data/Output" #Name directory where
10 studyname="Data" #name of study, only needed if datadir is a list of filenames
11 #=====
12
13 library(GGIR)
14
15 g.shell.GGIR(#-----
16             # General parameters
17             #-----
18             mode=mode, #specify above
19             datadir=datadir, #specify above
20             outputdir=outputdir, #specify above
21             studyname=studyname, #specify above
22             f0=f0, #specify above
23             f1=f1, #specify above
24             overwrite = TRUE, #overwrite previous milestone data?
25             do.imp=TRUE, # Do imputation? (recommended)
26             idloc=2, #id location (1 = file header, 2 = filename)
27             print.filename=TRUE,
28             storefolderstructure = FALSE,
29             #-----
30             # Part 1 parameters:
31             #-----
32             # Key functions: reading file, auto-calibration, and extracting features
33             window sizes = c(5,900,3600), #Epoch length, non-wear detection resolution, non-wear detection evaluation window
34             do.cal= TRUE, # Apply autocalibration? (recommended)
35             do.enmo = TRUE, #Needed for physical activity analysis
36             do.anglez=TRUE, #Needed for sleep detection
37             chunksize=1, #size of data chunks to be read (value = 1 is maximum)
38             printsummary=TRUE,
39             #-----
40             # Part 2 parameters:
41             #-----
42             # Key functions: Non-wear detection, imputation, and basic descriptives
43             strategy = 2, #Strategy (see tutorial for explanation)
44             ndayswindow=7, #only relevant when strategy = 3
45             hrs.del.start = 0, # Only relevant when strategy = 2. How many HOURS need to be ignored at the START of the measurement?
46             hrs.del.end = 0, # Only relevant when strategy = 2. How many HOURS need to be ignored at the END of the measurement?
47             maxdur = 8, # How many DAYS of measurement do you maximally expect?
48             includedaycrit = 16, # number of minimum valid hours in a day to attempt physical activity analysis
49             MSL5res = 10, #resolution in minutes of M5 and L5 calculation
50             winhr = c(5, 10), # size of M5 and L5 (5 hours by default)
51             qlevels = c(c(1380/1440),c(1410/1440)), #quantiles to calculate, set value at c() if you do not want quantiles
52             qwindow=c(0,8,21,24), #window over which to calculate quantiles
53             ilevels = c(), #c(0,100,400,8000), #acceleration values (metric ENMO) from which a frequency distribution needs to be
54             iglevels = TRUE,
55             mvpathreshold =c(100), #MVPA (moderate and vigorous physical activity threshold)
56             bout.metric = 4,
57             #-----
58             # Part 3 parameters:
59             #-----
60             # Key functions: Sleep detection
61             anglethreshold=5,
62             timethreshold=5,
63             ignorenonwear = TRUE, # if TRUE non-wear is not detected as sleep (if FALSE then it will work with imputed data)
64             do.part3.pdf = TRUE,
65             )
cc
24:36 (Untitled) R Script
Environment History Connections
```

```
mainscript_GGIR_demo.R* x
Source on Save Run Source
121 qwindow=c(0,8,21,24), #window over which to calculate quantiles
122 ilevels = c(), #c(0,100,400,8000), #acceleration values (metric ENMO) from which a frequency distribution needs to be
123 iglevels = TRUE,
124 mvpathreshold =c(100), #MVPA (moderate and vigorous physical activity threshold
125 bout.metric = 4,
126 #-----
127 # Part 3 parameters:
128 #-----
129 # Key functions: sleep detection
130 anglethreshold=5,
131 timethreshold= 5,
132 ignorenonwear = TRUE, # if TRUE non-wear is not detected as sleep (if FALSE then it will work with imputed data)
133 do.part3.pdf = TRUE,
134 #-----
135 # Part 4 parameters:
136 #-----
137 # Key functions: Integrating sleep log (if available) with sleep detection, storing day and person specific summaries
138 excludefirstlast = FALSE, # Exclude first and last night for sleep analysis?
139 includenightcrit = 16, # number of minimum valid hours in a day to attempt sleep analysis
140 def.noc.sleep = c(1), #c(),
141 # If sleep log is available:
142 # loglocation= "C:/Users/FlorianHuber/OneDrive - Netherlands eScience Center/Sprint_project_GGIR/Test-Data/sleeplogw.
143 outliers.only = FALSE,
144 criterror = 4,
145 relyonsleeplog = FALSE,
146 sleeplogidnum = TRUE, # Is the participant in the sleep log stored as a number (TRUE) or as a character (FALSE)
147 colid=1, #column in which the participant id or filename is stored
148 coln1=2, #column number for first day
149 do.visual = TRUE,
150 nnights = 9, #number of nights in the sleep log
151 #-----
152 # Part 5 parameters:
153 #-----
154 # Key functions: Merging physical activity with sleep analyses
155 excludefirstlast.part5 = FALSE,
156 save_ms5rawlevels = FALSE,
157 threshold.lig = c(30), #40 #threshold(s) for inactivity (can be more than one number)
158 threshold.mod = c(100), #100 120 #threshold(s) for moderate activity (can be more than one number)
159 threshold.vig = c(400), #500 #threshold(s) for vigorous activity (can be more than one number)
160 boutcriter = 0.8,
161 boutcriter.in = 0.9, #0.8 #fraction of an inactivity bout that needs to be below the threshold (needs to be 1 number)
162 boutcriter.lig = 0.8, #fraction of an light activity bout that needs to be between the thresholds (needs to be 1 numb
163 boutcriter.mvpa = 0.8, #fraction of an light activity bout that needs to be above the threshold (needs to be 1 number
164 boutdur.in = c(1,10,30), # duration of bouts to be calculated
165 boutdur.lig = c(1,10), # duration of bouts to be calculated
166 boutdur.mvpa = c(10), # duration of bouts to be calculated
167 timewindow = c("ww"), #,
168 #-----
169 # Report generation
170 #-----
171 # Key functions: Generating reports based on meta-data
172 do.report=c(2,4,5), #for what parts does and report need to be generated? (option: 2, 4 and 5)
173 visualreport=TRUE,
174 dofirstpage = TRUE, #first page of pdf-report with simple summary histograms
175 viewingwindow=2) #viewingwindow of visual report: 1 centres at day and 2 centers at night
```


Thank you and d**R**ink up!

Notes

- <https://www.youtube.com/watch?v=RuFBCAqFJ2M>
- <https://www.youtube.com/watch?v=S8YPTrYNWdU&feature=youtu.be>